

# Cicatrizial

Medical device that works by protecting the wound that has damaged the dermis, creating an optimal environment for the repair processes and at the same time protecting the area from external physical and bacterial agents.

MEDICAL DEVICE

CLASS II B

BIBLIOGRAPHIC  
EVIDENCE AVAILABLE

CICATRIZIAL performs an adjuvant action in the tissue repair process in the case of:

- Irritation and redness (including from diapers);
- Post peeling, hair removal and laser interventions;
- Navel of the newborn;
- Superficial injuries: cracks, scratches, grazes, abrasions, minor burns, superficial cuts;
- Deep wounds: surgical wounds, bedsores, ulcers

## DOSAGE AND INSTRUCTIONS FOR USE

Wash your hands before use. After cleaning and disinfecting the area to be treated, apply a thin layer of gel directly to the area to be treated.

In case of deep wounds, perform a light massage and apply a gauze on the injured part with a sterile bandage.

Apply twice a day, until symptoms disappear. The treatment should not exceed 30 days.

Dermatology



PACKAGING

25 g

GEL

Retail public price

in Italy: € 14.00

**DERMATOLOGICALLY TESTED  
FOR WOUNDS, BURNS AND IRRITATIONS**

## INGREDIENTS

ALOE VERA, SODIUM HYALURONATE, ACRYLATES/C10 – 30 ALKYL ACRYLATE CROSSPOLYMER, JOJOBA OIL\*, POLYVINYLPIRROLIDONE (PVP), SODIUM HYDRATE PEARLS, BENZYL ALCOHOL, GLYCERYL CAPRYLATE, GLYCERYL UNDECYLENATE, WATER

\* from organic farming